2101-C Nobel Ste. Julie, QC Canada, J3E 1Z8 p: 1.866.496.4777 f: 1.450.796.3888 salus-scientific.com

ActeroTM Salmonella/STEC Enrichment Media Product Information

INTENDED USE:

ActeroTM Salmonella/STEC Enrichment Media is a selective medium optimized for an improved enrichment of *Salmonella* spp. from food and environmental surface samples. Furthermore, the medium allows the enrichment of the Shiga toxin-producing *E. voli* (STEC) along or coenriched with *Salmonella* in ground beef.

PRINCIPLE OF OPERATION:

The principle of ActeroTM Salmonella/STEC Enrichment Media is based on the ability of *Salmonella* strains to optimize growth by the use of specific nutrients that are contained within the ActeroTM media. This unique medium formulation confers an important growth advantage when other bacteria are present.

KIT CONTENTS:

The kit contains sufficient material to prepare 35 liters of liquid medium.

- ✓ Dehydrated ActeroTM Salmonella/STEC Enrichment Media, bottle of 500 g.
- ✓ ActeroTM Salmonella/STEC Supplement # 1, 2 bottles of 35 mL.
- ✓ Actero™ Salmonella/STEC Supplement # 2, 1 bottle of 17 mL.

ADDITIONAL MATERIALS REQUIRED:

- 1. Distilled/deionized, sterile water.
- Sterile filter and non-filter Stomacher® bags.
- 3. Serological pipette, sterile.
- 4. Water bath 39-40 \pm 0.5°C
- 5. Incubator: at 39 \pm 0.5°C, 35 \pm 2°C
- 6. Tips and Adjustable Volume Pipette (100 1000 μL).
- 7. 10 µL calibrated inoculating loop
- 8. Rappaport-Vassiliadis Broth (RV).
- 9. Tetrathionate broth.
- 10. Xylose Lysine Tergitol-4 Agar (XLT4)
- 11. BG Sulfa Agar (BGS)
- 12. Xylose Lysine Deoxycholate agar (XLD)
- 13. Hektoen Enteric agar (HE)
- 14. Rapid'Salmonella
- 15. CHROMagarTM
- 16. Nonfat Dry Milk (NFDM)
- 17. Malachite green
- Stomacher® 3500/Stomacher® 400 (optional) available from multiple sources or equivalent.
- 19. Other regular laboratory equipment could also be required.

Environmental Samples

Non-bactericidal sterile cellulose sampling sponges (8×4×0.3 cm) premoistened with neutralizing Dey-Engley buffer (D/E). (Salus Scientific Inc. Cat # FCLS-005).

PROCEDURE:

ActeroTM Salmonella/STEC Enrichment Media Preparation

With the use of AUTOCLAVE

- Always shake the 500 g dry powder medium container before each use.
- 2. Measure 14.2 g of dry medium powder on the weight scale.

- 3. Suspend and mix this 14.2 g of the medium into a clean one-liter bottle of distilled water.
- 4. Sterilize it by autoclaving at 121°C for 15 min.
- Cool to room temperature and store at room temperature or refrigerate until use.
- Prior to use, the medium must be pre-warmed followed by the addition of 2 mL of supplement number 1 and 500 μL of supplement number 2 per liter of the medium.
- . The pH should be at 8.2 ± 0.2 prior to use.

Without the use of AUTOCLAVE

- 1. Always shake the 500 g dry powder medium container before each
- 2. Measure 14.2 g of dry medium powder on the weight scale.
- Suspend and mix this 14.2 g in one liter of pre-warmed <u>sterile</u> distilled water.
- The medium equilibrated now requires the addition of 2 mL of supplement number 1 and 0.5 mL of supplement number 2 per liter of the medium.
- 5. The medium prepared should be used immediately.
- 6. The pH should be at 8.2 ± 0.2 prior to use.

ActeroTM Salmonella/STEC Enrichment with 5% Nonfat Dry Milk (NFDM) Media Preparation

- 1. Always shake the 500 g dry powder medium container before use.
- 2. Measure 14.2 g of dry medium powder on the weight scale.
- 3. Suspend and mix this 14.2 g of the medium into a clean one-liter bottle of distilled water.
- 4. Add **50 g** of NFDM to the medium mixture.
- 5. Sterilize this bottle by autoclaving at 121°C for <u>6 min.</u>
- Remove the broth from the autoclave as soon as possible and cool to room temperature using ice water bath. Store refrigerated until use.
- Prior to use, the medium must be pre-warmed followed by the addition of 2 mL of supplement number 1 and 500 μL of supplement number 2 per liter of the medium.
- 8. The pH should be at 8.2 ± 0.2 prior to use.

Environmental Surface Sample Preparation

ActeroTM Salmonella Culture Method

- 1. Add to the non-bactericidal, non-bacteriostatic 8×4×0.3 cm sterile cellulose sampling sponge pre-moistened with D/E.
- 2. Wipe the surface to be tested with one side of the sponge (with excess liquid gently squeezed out) in a horizontal direction (approximately 10 cm), and with the other side in a vertical direction (approximately 10 cm) back and forth (one stroke back and one stroke forward) to cover the entire area of 100 cm².
- 3. Place each surface sampled sponge in a sterile sample bag, and keep it at 4 ± 2 °C until it is ready for testing. Analyze sample units as soon as possible after their reception in the laboratory.
- 4. When ready to test, pre-warm the prepared ActeroTM Salmonella Enrichment Media at $35 \pm 2^{\circ}C$ or $39 \pm 0.5^{\circ}C$.
- Add 90 ± 5 mL of the pre-warmed ActeroTM Salmonella Enrichment Media to each sponge sample in its sample bag.
- 6. Homogenize the sample for 30 seconds in a Stomacher® 400 or equivalent. Hand mixing is an acceptable alternative for stomaching. To hand mix, massage each sponge that is in the sealed Stomacher® bag for approximately 1 minute.



Summary of Enrichment and Analysis of Samples Using Actero™ Salmonella/STEC Enrichment Media

Sample Type	Sample Preparation	Analysis of Enriched Samples
Food Contact	 Swab a 100 cm² surface with a sponge (pre-moistened with D/E buffer) and keep it in a sterile bag at 4°C until tested. Homogenize the sample with 90 mL pre-warmed (39°C) ActeroTM Salmonella Enrichment Media. 	✓ By Actero™ Salmonella Culture Method
Surface > Stainless steel > Plastic	 Incubate at 39 ± 0.5°C for 18 h. Swab a 100 cm² surface with a sponge (pre-moistened with D/E buffer) and keep it in a sterile bag at 4°C until tested. Homogenize the sample with 90 mL pre-warmed (35°C) ActeroTM Salmonella Enrichment Media. Incubate at 35 ± 2°C for 14 - 18 h. 	✓ By BAX® System Real-Time PCR Assay for Salmonella ✓ By Actero™ Salmonella Culture Method
Non-Food Contact Surface Rubber Sealed concrete Ceramic	 Swab a 100 cm² surface with a sponge (pre-moistened with D/E buffer) and keep it in a sterile bag at 4°C until tested. Homogenize the sample with 90 mL pre-warmed (39°C) ActeroTM Salmonella Enrichment Media. Incubate at 39 ± 0.5°C for 18 h. 	✓ By Actero [™] Salmonella Culture Method
Ground Chicken	 Homogenize (30 s) 25 g sample with 50 mL pre-warmed (39°C) ActeroTM Salmonella Enrichment Media. Incubate at 39 ± 0.5°C for 20 h. 	✓ By Actero [™] Salmonella Culture Method
	 Homogenize (30 s) 25 g sample with 225 mL pre-warmed (35°C) ActeroTM Salmonella Enrichment Media supplemented with 50 mg/L malachite green. Incubate at 35 ± 2.0°C for 14-18 h. 	 ✓ By BAX® System Real-Time PCR Assay for Salmonella ✓ By Actero™ Salmonella Culture Method
Chicken Carcass Rinse	 Homogenize (30 s) 30 mL sample with 30 mL pre-warmed (35°C) Actero™ Salmonella Enrichment Media supplemented with 20 mg/L malachite green. Incubate at 35 ± 2.0°C for 16-20 h. 	✓ By BAX® System Real-Time PCR Assay for Salmonella ✓ By Actero™ Salmonella Culture Method
Ground Beef	 Homogenize (30 s) 325 g sample with 650 mL pre-warmed (39°C) ActeroTM Salmonella Enrichment Media. Incubate at 39.5 ± 0.5°C for 7 h in a water bath Transfer 0.5 mL of enriched sample into 10 mL TBG and 0,1 mL into 10 mL RVS and incubate at 42 ± 0.5°C for 22-24 h in a water bath. 	✓ By Actero [™] Salmonella Culture Method
	 Homogenize (60 s) 25 g sample with 75 mL pre-warmed (39°C) Actero™ Salmonella Enrichment Media supplemented with 50 mg/L malachite green. Incubate at 39 - 42± 0.5°C for 16-20 h. 	 ✓ By BAX® System Real-Time PCR Assay for Salmonella ✓ By Actero™ Salmonella Culture Method
	 Homogenize (60 s) 375 g sample with 1125 mL pre-warmed (39°C) Actero™ Salmonella Enrichment Media supplemented with 25 mg/L malachite green. Adjust pH to 7.0 ± 0.2 if necessary. Incubate at 39 ± 0.5°C for 20-24 h. 	 ✓ By BAX® System Real-Time PCR Assay for Salmonella ✓ By Actero™ Salmonella Culture Method
	 Homogenize (30 s) 100 g sample with 300 mL pre-warmed (39°C) ActeroTM Salmonella Enrichment Media. Incubate at 39 ± 0.5°C for 7 h in a water bath. 	✓ By Actero [™] Salmonella Culture Method
Whole Liquid Egg	 Homogenize (30 s) 100 g ± 2 g sample with 700 mL pre-warmed (39°C) ActeroTM Salmonella Enrichment Media. Incubate at 39 ± 0.5°C for 18 h. 	✓ By Actero™ Salmonella Culture Method
	 Homogenize (30 s) 100 g ± 2 g sample with 300 mL pre-warmed (35°C) ActeroTM Salmonella Enrichment Media. Adjust pH to 7.0 ± 0.2 if necessary. Incubate at 35 ± 2.0°C for 18-22 h. 	 ✓ By BAX® System Real-Time PCR Assay for Salmonella ✓ By Actero™ Salmonella Culture Method
Dried Whole Egg	 Homogenize (30 s) 100 g sample with 300 mL pre-warmed (39°C) ActeroTM Salmonella Enrichment Media supplemented with 5% NFDM. Incubate at 35 ± 2.0°C for 14-18 h. 	✓ By Actero™ Salmonella Culture Method
1-88	 Homogenize (30 s) 100 g sample with 600 mL pre-warmed (35°C) Actero™ Salmonella Enrichment Media supplemented with 5% NFDM. Incubate at 35 ± 2.0°C for 14-18 h. 	✓ By BAX® System Real-Time PCR Assay for Salmonella ✓ By Actero™ Salmonella Culture Method
Shell Egg	 Homogenize by hand a sample of 20 eggs with 1000 mL pre-warmed (35°C) Actero™ Salmonella Enrichment Media. Incubate at 35 ± 2°C for 16-20 h. 	✓ By BAX® System Real-Time PCR Assay for Salmonella ✓ By Actero™ Salmonella Culture Method
Raw Frozen Scallop	 Homogenize (30 s) 25 g sample with 50 mL pre-warmed (39°C) Actero™ Salmonella Enrichment Media. For water bath: Incubate at 39 ± 0.5°C for 7 h. For incubator: Incubate at 39 ± 0.5°C for 18 h. 	✓ By Actero [™] Salmonella Culture Method
Sprout	 Homogenize (60 s) 25 g sample with 150 mL pre-warmed (39°C) Actero™ Salmonella Enrichment Media. Incubate at 39 ± 0.5°C for 7 h in a water bath. Transfer 1,0 mL of enriched sample into 10 mL TBG and 0,1 mL into 10 mL RVS and incubate respectively at 43 ± 0.2°C C and at 42 ± 0.2°C for 18 h in a water bath. 	✓ By Actero™ Salmonella Culture Method

2101-C Nobel Ste. Julie, QC Canada, J3E 1Z8 p: 1.866.496.4777 f: 1.450.796.3888 salus-scientific.com

Sample Type	1. Sample Preparation	✓ Analysis of Enriched Samples
Milk Chocolate	2. Homogenize (60 s) 25 g sample with 175 mL pre-warmed (35°C) Actero™ Salmonella	✓ By BAX® System Real-Time PCR Assay
	Enrichment Media. 3. Incubate at 35 ± 2°C for 22-26 h.	for Salmonella ✓ By Actero™ Salmonella Culture Method
		✓ By Actero [™] Salmonella Culture Method
	 Homogenize (60 s) 25 g sample with 175 mL pre-warmed (35°C) Actero™ Salmonella Enrichment Media. 	✓ By Actero™ Salmonella Culture Method
	2. Incubate at 39 ± 0.5°C for 18 h.	by Actero Samonena Culture Method
Chocolate Liquor	Homogenize (60 s) 25 g sample with 225 mL pre-warmed (35°C) Actero™ Salmonella	✓ By BAX® System Real-Time PCR Assay
	Enrichment Media.	for Salmonella
	2. Incubate at $35 \pm 2^{\circ}$ C for 26-30 h.	✓ By Actero™ Salmonella Culture Method
Cocoa Powder	1. Homogenize (60 s) 25 g sample with 175 mL pre-warmed (35°C) Actero™ Salmonella	✓ By BAX® System Real-Time PCR Assay
	Enrichment Media supplemented with 5% NFDM.	for Salmonella
	2. Incubate at 35 ± 2°C for 16-20 h.	✓ By Actero [™] Salmonella Culture Method
	1. Homogenize (60 s) 25 g sample with 225 mL pre-warmed (35°C) Actero TM Salmonella	✓ By BAX® System Real-Time PCR Assay
	Enrichment Media.	for Salmonella
Dry Pet Food	2. Incubate at $35 \pm 2^{\circ}$ C for $18-22$ h.	✓ By Actero [™] Salmonella Culture Method
Diy 1 ct 1 ood	1. Homogenize (60 s) 375 g sample with 2625 mL pre-warmed (35°C) Actero TM Salmonella	✓ By BAX® System Real-Time PCR Assay
	Enrichment Media.	for Salmonella
	2. Incubate at 35 ± 2°C for 18-22 h.	✓ By Actero™ Salmonella Culture Method
D 41 1	1. Homogenize (60 s) 375 g sample with 750 mL pre-warmed (35°C) Actero™ Salmonella	✓ By BAX® System Real-Time PCR Assay
Raw Almond	Enrichment Media.	for Salmonella
	2. Incubate at 35 ± 2°C for 16-20 h.	✓ By Actero™ Salmonella Culture Method
Peanut Butter	1. Homogenize (60 s) 25 g sample with 175 mL pre-warmed (35°C) Actero™ Salmonella Enrichment Media.	✓ By BAX® System Real-Time PCR Assay for Salmonella
	2. Incubate at 35 ± 2°C for 16-20 h.	✓ By Actero™ Salmonella Culture Method
Dried Parsley	Homogenize (60 s) 25 g sample with 225 mL pre-warmed (35°C) Actero TM Salmonella	✓ By BAX® System Real-Time PCR Assay
	Enrichment Media then manually to remove parsley from the bag walls.	for Salmonella
	2. Incubate at 35 ± 2°C for 20-24 h.	✓ By Actero™ Salmonella Culture Method
Dried Raisin	Homogenize (60 s) 25 g sample with 75 mL pre-warmed (35°C) Actero TM Salmonella	✓ By BAX® System Real-Time PCR Assay
	Enrichment Media.	for Salmonella
	2. Incubate at 35 ± 2°C for 16-20 h.	✓ By Actero™ Salmonella Culture Method
Whole Black Pepper	1. Homogenize 25 ± 1 g sample manually with 75 ± 5 mL of pre-warmed Actero™ Salmonella	✓ By BAX® System Real-Time PCR Assay
	Enrichment Media.	for Salmonella
	2. Incubate at 35 ± 2°C for 16 -20 h.	✓ By Actero [™] Salmonella Culture Method

Enrichment of Environmental Surface Samples Using ActeroTM Salmonella/STEC Enrichment Media

Stainless Steel, Plastic

ActeroTM Salmonella Culture Method

For the enrichment phase, close the bag and incubate the sample in an incubator for 18 ± 0.5 h at 39 ± 0.5 °C. Adherence to temperature is important for accurate results.

$\underline{BAX \circledast}$ System Real-Time PCR Assays for Salmonella or Actero Salmonella Culture Method

For the enrichment phase, close the bag loosely and incubate the sample in an incubator for $16 \pm 2 h$ at $35 \pm 2^{\circ}C$.

Rubber, Ceramic Tile and Sealed Concrete

ActeroTM Salmonella Culture Method

For the enrichment phase, close the bag and incubate the sample in an incubator for 18 ± 0.5 h at 39 ± 0.5 °C. Adherence to temperature is important for accurate results.

At the end of the enrichment period, mix the sample thoroughly and transfer 10.0 ± 0.1 mL of the enriched sample to a tube. Cap the tube.

Preparation and Enrichment of Food Samples Using ActeroTM Salmonella/STEC Enrichment Media

Raw Ground Chicken (25 g)

ActeroTM Salmonella Culture Method

 Add 50 mL of pre-warmed (39°C) and supplemented ActeroTM Salmonella broth to 25 g of sample in a filter-equipped Stomacher® bag.

- Homogenize the sample for 30 seconds in a Stomacher[®] 400 circulator or equivalent or mix vigorously in the Stomacher[®] bag for 1 minute if there is no Stomacher[®] machine available.
- Close the bag loosely and incubate the sample upright for 20 h at 39°C in an incubator for enrichment.
- After 20 hours, remove the sample from the incubator, re-suspend the contents by shaking the bag and transfer 10.0 ± 0.1 mL to a tube. Cap the tube.

BAX® System Real-Time PCR Assay for Salmonella or ActeroTM Salmonella Culture Method

- 1. Pre-warm the ActeroTM Salmonella broth supplemented with 50 mg/L malachite green at $35 \pm 2^{\circ}\text{C}$ before use.
- Add 225 ± 5 mL of pre-warmed ActeroTM Salmonella broth to each filter bag containing the 25 ± 1 g test portion.
- 3. Mix each sample by hand to homogenize it.
- Close the bag loosely, and incubate the sample for 14 18 h at 35 ± 2°C using an incubator.

Chicken Carcass Rinse

BAX® System Real-Time PCR Assays for Salmonella or ActeroTM Salmonella Culture Method

- 1. Pre-warm the ActeroTM Salmonella broth supplemented with 20 mg/L malachite green at $35 \pm 2^{\circ}C$ before use.
- 2. Add 30 ± 1 mL of pre-warmed ActeroTM Salmonella broth to each filter bag containing the 30 ± 1 mL test portion.
- 3. Mix each sample by hand to homogenize it.
- Close the bag loosely, and incubate the sample for 16 − 20 h at 35 ± 2°C using an incubator.



2101-C Nobel Ste. Julie, QC Canada, J3E 1Z8 p: 1.866.496.4777 f: 1.450.796.3888

salus-scientific.com

Raw Ground Beef (325 g)

ActeroTM Salmonella Culture Method

- Add 650 mL of pre-warmed and supplemented Actero TM Salmonella broth to 325 g of sample in a filter-equipped Stomacher® bag.
- Homogenize the sample for 30 seconds in a Stomacher® 3500 or equivalent. Alternatively, mix vigorously in the bag for 1 minute if there is no Stomacher® machine available.
- 3. Close the bag loosely and incubate the samples for 7 h at 39.5°C in a water bath for enrichment. If there are a large number of samples to be analyzed, verify that the temperature of the water between the sample bags reaches 39.5°C before starting to record the required incubation time. It is important to precisely control the enrichment period to obtain valuable accurate results.
- After 7 hours, remove the sample from the water bath, mix the contents by shaking the bag and transfer 10.0 ± 0.1 mL to a tube. Cap the tube.
- Transfer 0.5 ml of enriched sample into 10 mL Tetrathionate Broth and 0.1 ml into 10 mL modified Rappaport-Vassiliadis broth and incubate tubes at 42 ± 0.5 °C for 22-24 h.

Note: Test limitation: test standardization is for meat with a maximum aerobic microbial flora of 4 x 105 CFU/g

Raw Ground Beef (25 g)

BAX® System Real-Time PCR Assay for Salmonella or ActeroTM Salmonella Culture Method

- 1. Pre-warm the ActeroTM Salmonella broth supplemented with 50 mg/L malachite green at 39-42 \pm 0.5°C before use.
- Add $75 \pm 5 \text{ mL}$ of pre-warmed ActeroTM Salmonella broth to each filter bag containing the $25 \pm 1\,g$ test portion.
- Homogenize each sample for 60 s using a stomacher.
- Close the bag loosely, and incubate the 25 g sample for 16 20 h at $39-42 \pm 0.5$ °C using an incubator.

Raw Ground Beef (375 g)

BAX® System Real-Time PCR Assay for Salmonella or ActeroTM Salmonella Culture Method

- 1. Pre-warm the ActeroTM Salmonella broth supplemented with 25 mg/L malachite green at 39 ± 0.5 °C before use.
- Add 1125 ± 25 mL of pre-warmed Actero TM Salmonella broth to each filter bag containing the 375 \pm 10 g test portion.
- Homogenize each sample for 60 s using a stomacher.
- Adjust pH to 7.0 ± 0.2 .
- Close the bag loosely, and incubate the 375 g for 20 24 h at 39 ± 0.5 °C using an incubator.

Whole Liquid Egg (100 g)

ActeroTM Salmonella Culture Method: 7 h Enrichment in Water Bath

- 1. Add 300 ml of pre-warmed (39°C) and supplemented ActeroTM Salmonella broth to 100 g of sample in a filter-equipped Stomacher® bag. Adjust pH, if necessary, to 7.0 ± 0.4 .
- Homogenize the sample for 30 seconds at 150 rpm in a Stomacher® 3500 or equivalent. Alternatively, mix in the bag until homogeneous if there is no Stomacher® machine available.
- 3. Close the bag loosely and incubate the sample upright for 7 h at 39°C in a water bath for enrichment. If there are a large number of samples to be analyzed, verify that the temperature of the water between the sample bags reaches 39°C before starting to record the required incubation time. It is important to precisely control the enrichment period to obtain valuable accurate results.
- After 7 hours, remove the sample from the water bath, mix the contents by shaking the bag and transfer $10.0 \pm 0.1 \text{ mL}$ to a tube. Cap the tube.

ActeroTM Salmonella Culture Method: 18 h Enrichment in Incubator

- Add 700 ml of pre-warmed (39°C) and supplemented ActeroTM Salmonella broth to 100 g of sample in a non-filtered Stomacher® bag. Adjust pH, if necessary, to 7.0 ± 0.4 .
- Homogenize the sample for 30 seconds in a Stomacher® 3500 or equivalent. Alternatively, mix vigorously in the bag for 1 minute if there is no Stomacher® machine available.
- Close the bag loosely and incubate the sample upright for 18 h at 39°C in an incubator for enrichment.
- After 18 hours, remove the sample from the incubator, mix the contents by shaking the bag and transfer 10.0 ± 0.1 mL to a tube. Cap the tube.

BAX® System Real-Time PCR Assay for Salmonella or ActeroTM Salmonella Culture Method

- Add 300 ml of pre-warmed (35°C) and supplemented ActeroTM Salmonella broth to 100 g of sample in a non-filtered Stomacher®
- Mix vigorously by hand in the bag to homogenize it.
- 3. Adjust pH to 7.0 ± 0.2 .
- Close the bag loosely and incubate the sample upright for 18-22 h at 35°C in an incubator for enrichment.

Dried Whole Egg (100 g)

BAX® System Real-Time PCR Assay for Salmonella or ActeroTM Salmonella Culture Method

- Pre-warm the ActeroTM Salmonella broth supplemented with 5% NFDM at $35 \pm 2^{\circ}$ C before use.
- Add 600 ± 10 mL of pre-warmed ActeroTM Salmonella broth supplemented with 5% NFDM to each filter bag containing the 100 \pm 5 g test portion.
- 3. Mix each sample by hand to homogenize it.
- 4. Close the bag loosely, and incubate the sample for 14 18 h at 35 \pm 2°C using an incubator.

Shell Egg

BAX® System Real-Time PCR Assay for Salmonella or ActeroTM Salmonella Culture Method

- Add $1000 \pm 50 \text{ mL}$ of pre-warmed $(35 \pm 2^{\circ}\text{C})$ ActeroTM Salmonella broth to each filter bag containing the 20-egg test portion.
- Homogenize each sample by hand mixing.
- Close the bag loosely, and incubate the sample for 16-20 h at 35 \pm 2°C using an incubator

Raw Frozen Scallop (25 g)

ActeroTM Salmonella Culture Method

- Add 50 ml of pre-warmed (39°C) and supplemented ActeroTM Salmonella broth to 25 g of sample in a filter-equipped Stomacher®
- Homogenize the sample for 30 seconds in a Stomacher® 400 circulator or equivalent. Alternatively, mix vigorously in the bag for 1 minute if there is no Stomacher® machine available.

For a 7 h Enrichment in Water Bath

- Close the bag loosely and incubate the sample upright for 7 h at 39°C in a water bath for enrichment. If there are a large number of samples to be analyzed, verify that the temperature between the sample bags reaches 39°C before starting to record the incubation time. It is important to precisely control the enrichment period to obtain valuable and accurate results.
- After 7 hours, remove the samples from the water bath, mix the contents by shaking the bag and transfer 10.0 ± 0.1 mL to a tube. Cap

For a 18 h Enrichment in Incubator

Close bag loosely and incubate sample upright for 18 h at 39°C in an incubator for enrichment.



2101-C Nobel Ste. Julie, QC Canada, J3E 1Z8 p: 1.866.496.4777 f: 1.450.796.3888 salus-scientific.com

- 2. Homogenize each sample for **1 min** in a Stomacher® 400 circulator or equivalent.
- Close the bag loosely, and incubate the sample for 16-20 h at 35 ± 2°C using an incubator.

Sprout (25 g)

the tube.

ActeroTM Salmonella Culture Method

 Add 150 ml of pre-warmed (39°C) ActeroTM Salmonella broth to 25 g of sample in a filter-equipped Stomacher® bag.

After 18 hours, remove the samples from the incubator, mix the

contents by shaking the bag and transfer 10.0 ± 0.1 mL to a tube. Cap

- 2. Homogenize the sample for **60 seconds** in a Stomacher[®] 400 circulator or equivalent. Alternatively, mix vigorously in the bag for 1 minute if there is no Stomacher® machine available.
- 3. Close the bag loosely and incubate the sample upright for 7 h at 39°C in a water bath for enrichment. If a large number of samples are to be analyzed, verify that the temperature between the sample bags reaches 39°C before starting to record the incubation time. It is important to precisely control the enrichment period to obtain valuable and accurate results.
- 4. After 7 hours, remove the sample from the water bath, mix the contents by shaking the bag and transfer 10.0 ± 0.1 mL to a tube. Cap the tube
- 5. Transfer 1.0 ml of enriched sample into 10 mL Tetrathionate Broth and 0.1 ml into 10 mL Rappaport-Vassiliadis broth (RV) and incubate respectively at 43 ± 0.2°C and 42 ± 0.2°C for 18 h (because sprouts are considered to have high microbial load).

Dry Pet Food

BAX® System Real-Time PCR Assay for Salmonella or ActeroTM Salmonella Culture Method

- Add 225 ± 5 mL of pre-warmed (35 ± 2°C) Actero™ Salmonella broth to 25 g of sample in a filter-equipped Stomacher® bag.
- 2. Add 2625 ± 25 mL of pre-warmed Actero™ Salmonella broth to 375 g of sample in a filter-equipped Stomacher® bag.
- 3. Homogenize each sample for **60 seconds** in a Stomacher® 400 circulator or Stomacher® 3500 or equivalent.
- Close the bag loosely, and incubate the sample for 18-22 h at 35 ± 2°C using an incubator.

Milk Chocolate

- 1. Pre-warm the ActeroTM Salmonella medium before use.
- 2. Add 175 ± 5 mL of pre-warmed (35 ± 2°C) ActeroTM Salmonella broth to 25 g of sample in a filter-equipped Stomacher® bag.
- Homogenize each sample for 60 seconds in a Stomacher® 400 circulator or equivalent.

BAX® System Real-Time PCR Assay for Salmonella

 Close the bag loosely, and incubate the sample for 22-26 h at 35 ± 2°C using an incubator.

ActeroTM Salmonella Culture Method

5. Close the bag loosely, and incubate the sample for 18-22 h at 39 ± 0.5°C or 22-26 h at 35 ± 2°C using an incubator.

Chocolate Liquor

BAX® System Real-Time PCR Assay for Salmonella or ActeroTM Salmonella Culture Method

- Add 225 ± 5 mL of pre-warmed (35 ± 2°C) Actero™ Salmonella broth to 25 g of sample in a filter-equipped Stomacher® bag.
- Homogenize each sample for 2 minutes at 175 rpm in a Stomacher® 400 circulator or equivalent.
- 3. Close the bag loosely, and incubate the sample for 26-30 h at 35 \pm 2°C using an incubator.

Cocoa Powder

BAX® System Real-Time PCR Assay for Salmonella or ActeroTM Salmonella Culture Method.

Add 175 ± 5 mL of pre-warmed (35 ± 2°C) Actero™ Salmonella broth supplemented with 5% NFDM to each filter bag containing the 25 ± 1 g test portion.

Raw Almond

BAX® System Real-Time PCR Assay for Salmonella or ActeroTM Salmonella Culture Method

- 1. Pre-warm the ActeroTM Salmonella broth at $35 \pm 2^{\circ}$ C before use.
- Add 750 ± 25 mL of pre-warmed Actero™ Salmonella broth to each filter bag containing the 375 ± 10 g sample.
- 3. Homogenize each sample for **60 s** using a stomacher.
- Close the bag loosely, and incubate the sample for 16 20 h at 35 ± 2°C using an incubator.

Peanut Butter

BAX® System Real-Time PCR Assay for Salmonella or ActeroTM Salmonella Culture Method

- 1. Pre-warm the ActeroTM Salmonella broth at $35 \pm 2^{\circ}$ C before use.
- 2. Add 175 ± 5 mL of pre-warmed ActeroTM Salmonella broth to each filter bag containing the 25 ± 1 g test portion.
- 3. Homogenize each sample for 60 s using a stomacher.
- Close the bag loosely, and incubate the sample for 16 20 h at 35 ± 2°C using an incubator.

Dried Parsley

BAX® System Real-Time PCR Assay for Salmonella or ActeroTM Salmonella Culture Method

- 1. Pre-warm the ActeroTM Salmonella broth at $35 \pm 2^{\circ}$ C before use.
- 2. Add 225 ± 5 mL of pre-warmed ActeroTM Salmonella broth to each filter bag containing the 25 ± 1 g test portion.
- 3. Homogenize each sample for 60 s using a stomacher and then mix manually to remove the dried parsley from the sides of the bags.
- Close the bag loosely, and incubate the sample for 20-24 h at 35 ± 2°C using an incubator.

Dried Raisin

BAX® System Real-Time PCR Assay for Salmonella or ActeroTM Salmonella Culture Method

- 1. Pre-warm the ActeroTM Salmonella broth at $35 \pm 2^{\circ}$ C before use.
- 2. Add 75 ± 5 mL of pre-warmed ActeroTM Salmonella broth to each filter bag containing the 25 ± 1 g test portion.
- 3. Homogenize each sample for 60 s using a stomacher.
- Close the bag loosely, and incubate the sample for 16 -20 h at 35 ± 2°C using an incubator.

Whole Black Pepper

BAX® System Real-Time PCR Assay for Salmonella or ActeroTM Salmonella Culture Method

- 1. Pre-warm the Actero TM Salmonella broth at $\mathbf{35} \pm 2^{\circ}\mathbf{C}$ before use.
- 2. Add 75 \pm 5 mL of pre-warmed ActeroTM Salmonella broth to each filter bag containing the 25 \pm 1 g test portion.
- 3. Homogenize each sample manually.
- Close the bag loosely, and incubate the sample for 16 -20 h at 35 ± 2°C using an incubator.

Analysis of Enriched Samples

ActeroTM Salmonella Culture Method

Raw ground chicken (39°C), raw ground beef (39.5°), dried whole eggs, whole liquid eggs:

Streak the samples directly onto selective agar plates (XLT4 and BGS) using a calibrated loop of 10 μL and, if necessary, follow the confirmation procedure as recommended in the current USDA FSIS Microbiology Laboratory Guidebook chapter for Salmonella.



2101-C Nobel Ste. Julie, QC Canada, J3E 1Z8 p: 1.866.496.4777 f: 1.450.796.3888

salus-scientific.com

Raw ground chicken (35°C), chicken carcass rinse

Streak the enriched sample directly onto Rapid'Salmonella and BGS agar plates using a calibrated loop of 10 µL and, if necessary, follow the confirmation procedure as recommended in the current USDA FSIS Microbiology Laboratory Guidebook chapter for Salmonella.

Raw ground beef (25 g)

Streak the enriched sample directly onto CHROMAgar TM Salmonella and BGS agar plates using a calibrated loop of 10 µL and, if necessary, follow the confirmation procedure as recommended in the current USDA FSIS Microbiology Laboratory Guidebook chapter for Salmonella.

Raw ground beef (375 g)

✓ Streak the enriched sample directly onto Rapid'Salmonella and XLT4 agar plates using a calibrated loop of 10 µL and, if necessary, follow the confirmation procedure as recommended in the current USDA FSIS Microbiology Laboratory Guidebook chapter for Salmonella.

Raw frozen scallops, sprouts, dry pet food, milk chocolate, chocolate liquor, raw almond, peanut butter, dried raisin, dried parsley, whole black pepper, whole liquid egg, environmental surface samples:

Streak the samples onto selective agar plates (XLD and HE) using a calibrated loop of 10 µL and, if necessary, follow the confirmation procedure as recommended in the US FDA Bacteriological Analytical Manual Chapter 5.

Interpretation and Test Result Report

- ✓ All samples presenting typical colony (ies) after 24 or 48 h in the selective agar should be considered as presumptive positive(s). The presumptive results should be confirmed according to the US FDA Bacteriological Analytical Manual Chapter 5 and the current USDA FSIS Microbiology Laboratory Guidebook chapter for Salmonella.
- All samples which do not present typical colonies after 48h of incubation can be considered negative.

BAX® System Real-Time PCR Assay for Salmonella

Enriched Sample Preparation Prior to Testing

- 1. Label and arrange a 2 mL microcentrifuge tube in a rack.
- 2. Add 2 mL of PBS to each microcentrifuge tube and place it with an open cap in a microcentrifuge tube rack.
 - a. For raw ground beef, raw ground chicken, chicken carcass rinse, dried whole egg, whole liquid egg, dry pet food, milk chocolate, chocolate liquor, raw almond, peanut butter, dried raisin, dried parsley and whole black pepper, transfer 80 µL enriched sample from the bag into each microcentrifuge tube.
 - b. For cocoa powder, shell egg, stainless steel and plastic environmental sponge samples, transfer $40~\mu L$ enriched sample from the bag into each microcentrifuge tube.
- 3. Refer to the Test Protocol section of the package insert for the BAX® System Real-Time PCR Assay for Salmonella for the following steps.

PRODUCT STORAGE AND SHELF LIFE:

The dehydrated medium and supplement number 1 should be stored at room temperature (12-25°C), in a tightly closed bottle in a cool dry place. Supplement number 2 should be stored at refrigerated temperature (2-8°C) in a cool dry place protected from light. The expiration dates are indicated on the packaging.

The prepared autoclaved medium without supplement can be stored for up to 6 months and the **supplemented** medium can be stored for 2 months at 2-8°C and 1 month at 18-25°C, in a tightly closed bottle protected from light. Please take into consideration that the medium should be autoclaved and manipulated in aseptic conditions.

DISPOSAL:

Dispose of all materials used and the enrichment media by autoclaving or according to an approved practice. Ensure that all biohazardous waste is disposed of according to local, municipal, provincial, state and/or federal regulations.

PRECAUTIONS:

Salmonella are categorized as Biosafety Level 2 pathogens. Biosafety level 2 procedures should be exercised (https://www.cdc.gov/labs/BMBL.html). The use of microbiological media such as the Actero™ Salmonella Enrichment Media requires trained laboratory personnel familiar with good microbiological laboratory practices. Wear a laboratory coat, disposable gloves and eye protection while handling specimens and performing the assay is strongly recommended. Material Safety Data Sheet (MSDS) must be obtained from the manufacturer for the media, chemicals, reagents and microorganisms used in the analysis. The personnel who will handle the material should read the MSDS prior to start-up.

All enrichment broths may contain various pathogens whether they contain Salmonella spp. or not. Furthermore, some pathogen bacteria have a very low infective dose (Ex. E. coli O157:H7 is estimated to be 50 organisms). Thus, extreme care should be taken in handling test samples and enrichment broths.

TERMS AND CONDITIONS:

Salus Scientific Inc. makes no representations and warranties concerning its products other than those stated herein. All Product(s) delivered hereunder by Salus Scientific Inc., its affiliates or any other person on its behalf shall, at the time of delivery, be manufactured to meet Salus Scientific Inc.'s specifications and all applicable laws. All other terms, conditions and warranties, including any warranty of merchantability, quality, fitness or suitability for a particular or intended purpose, implied by common law or statue, (implied warranties) are expressly excluded.

CATALOGUE NUMBER:

FCM-010 (BBFCM-010): ActeroTM Salmonella/STEC Enrichment Media, 500 g

FOR FURTHER INFORMATION, PLEASE CONTACT:

Salus Scientific Inc. Suite 450, 1414 – 8 St. S.W. Calgary, Alberta, Canada T2R 1J6 Tel: 1-877-298-0208



